






Wkręty korowe samogwintujące Ø 1.5

Ø 1.1  REF

L(mm) REF

T4 6,7,8...16 **80-1-01-330-L**

16,18,20...30

WIERTŁO Ø1.1		80-40-23-11.93.AO
PROWADNIK WIERTŁA		80-44-06-15
RĘKOJEŚĆ		80-101-10
TRZPIEŃ		80-43-270-T4.70

Wkręty blokowane samogwintujące Ø 1.5

Ø 1.1  REF


L(mm) REF

T4 4,5,6...15 **80-1-01-329-L**

16,18,20...30

WIERTŁO Ø1.1		80-40-23-11.93.AO
PROWADNIK WIERTŁA		80-44-06-15
RĘKOJEŚĆ		80-101-10
TRZPIEŃ		80-43-270-T4.70





Wkręty blokowane zmiennokątowe Ø 1.5

Ø 1.1  REF

L(mm) REF

T4 4,5,6...16 **80-1-01-381-L**

16,18,20...30

WIERTŁO Ø1.1		80-40-23-11.93.AO
PROWADNIK WIERTŁA		80-44-07-15
RĘKOJEŚĆ		80-101-10
TRZPIEŃ		80-43-270-T4.70

Wkręty korowe Ø 2.4 samogwintujące



Ø 1.8  REF

L(mm) REF

T8 8,10,12...40 **80-X-01-240-L**

L(mm) REF

8,10,12...60 **80-X-01-239-L**

WIERTŁO Ø1.8		
PROWADNIK WIERTŁA		80-4-57-192-00
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		
		S2.5 80-28-5-3
	T8	80-43-268-T8.115

Wkręty blokowane Ø 2.4 samogwintujące

Ø 1.8  REF

L(mm) REF


T8 8,10,12...40 **80-X-01-208-L**

L(mm) REF

8,10,12...40 **80-X-01-200-L**


WIERTŁO Ø1.8		
PROWADNIK WIERTŁA		80-4-57-49-00
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		
		S2.5 80-28-5-3
	T8	80-43-268-T8.115

Wkręty blokowane zmiennokątowe Ø 2.4


Ø 1.8  REF

L(mm) REF

T8 8,10,12...40 **80-X-01-298-L**

WIERTŁO Ø1.8		
PROWADNIK WIERTŁA		80-4-57-171-01
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		80-43-268-T8.115

Wkręty korowe samogwintujące Ø 2.0

Ø 1.5  REF

L(mm) REF






T6 4,5,6...12 **80-X-01-282-L**

14,16,18...38


L(mm) REF

6,7,8...16 **80-X-01-45-L**

16,18,20...30

WIERTŁO Ø1.5		80-4-40-03-15
PROWADNIK WIERTŁA		80-4-57-191-00
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		
		S1.5 80-28-5-2
	T6	80-43-268-T6.115

Wkręty blokowane samogwintujące Ø 2.0

Ø 1.5  REF

L(mm) REF






T6 6,7,8...13 **80-X-01-262-L**

14,16,18...40


L(mm) REF

6,7,8...13 **80-X-01-260-L**

14,16,18...40

WIERTŁO Ø1.5		80-4-40-03-15
PROWADNIK WIERTŁA		80-57-122-00
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		
		S1.5 80-28-5-2
	T6	80-43-268-T6.115





Wkręty blokowane zmiennokątowe Ø 2.0

Ø 1.5  REF

L(mm) REF

T6 6,7,8...13 **80-X-01-380-L**

14,16,18...30

WIERTŁO Ø1.5		80-4-40-03-15
PROWADNIK WIERTŁA		80-4-57-212-00
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		80-43-268-T6.115

Wkręty korowe Ø 2.7 samogwintujące

Ø 2.0  REF

L(mm) REF

T8 8,10,12...60 **80-X-01-216-L**

L(mm) REF

8,10,12...60 **80-X-01-48-L**

WIERTŁO Ø2.0		
PROWADNIK WIERTŁA		80-4-57-192-00
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		
		S2.5 80-28-5-3
	T8	80-43-268-T8.115

Wkręty blokowane Ø 2.7 samogwintujące

Ø 2.0  REF

L(mm) REF


T8 10,12,14...60 **80-X-01-209-L**

L(mm) REF

8,10,12...60 **80-X-01-201-L**





WIERTŁO Ø2.0		
PROWADNIK WIERTŁA		80-4-57-74-00
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		
		S2.5 80-28-5-3
	T8	80-43-268-T8.115








































Wkręty blokowane zmiennokątowe Ø 2.7

Ø 2.0  REF

L(mm) REF

T8 10,16,18...40 **80-X-01-325-L**

WIERTŁO Ø2.0		
PROWADNIK WIERTŁA		80-4-57-171-02
RĘKOJEŚĆ		80-28-5-1
TRZPIEŃ		80-43-268-T8.115

Wkręty korowe samogwintujące  Ø 2.5 L(mm) REF		WIERŁO Ø2.5 		Wkręty blokowane samogwintujące  Ø 2.7 L(mm) REF		WIERŁO Ø2.7 	
T15 8,10,12...70 75,80,85...90 L(mm) REF	80-X-01-221-L	PROWADNIK WIERŁA 80-4-57-191-00 	T15 10,12,14...60 65,70,75...95 L(mm) REF	80-X-01-223-L	PROWADNIK WIERŁA 80-4-57-74-00 	RĘKOJEŚĆ 80-28-5-1 	RĘKOJEŚĆ 80-28-5-1 
T15 8,10,12...70 75,80,85...90 L(mm) REF	80-X-01-51-L	TRZPIEŃ 	T15 10,12,14...60 65,70,75...95 L(mm) REF	80-X-01-11-L	TRZPIEŃ 	 S2.5 80-28-5-3	 S2.5 80-28-5-3
		 T15 80-43-268-T15.115			 T15 80-43-268-T15.115		
Wkręty blokowane zmiennokątowe Ø 3.5  Ø 2.7 L(mm) REF		WIERŁO Ø2.7 		Wkręty gąbczaste samogwintujące  Ø 2.0 L(mm) REF		WIERŁO Ø2.0 	
T15 10,12,14...60 65,70,75...95 L(mm) REF	80-X-01-322-L	PROWADNIK WIERŁA 80-4-57-171-03 	T15 10,12,14...50 50,55,60...75 L(mm) REF	80-X-01-321-L	PROWADNIK WIERŁA 80-4-57-192-00 	RĘKOJEŚĆ 80-28-5-1 	RĘKOJEŚĆ 80-28-5-1 
		TRZPIEŃ 80-43-268-T15.115			TRZPIEŃ 	 S2.5 80-28-5-3	 S2.5 80-28-5-3
Wkręty korowe samogwintujące  Ø 3.2 L(mm) REF		WIERŁO Ø3.2 80-4-40-08-32 		Wkręty blokowane samogwintujące  Ø 3.7 L(mm) REF		WIERŁO Ø3.7 80-4-40-06-37 	
T25 12,14,16...90 90,95,100...115 L(mm) REF	80-X-01-218-L	PROWADNIK WIERŁA 80-4-57-57-00 	T25 14,16,18...60 65,70,75...95 L(mm) REF	80-X-01-222-L	PROWADNIK WIERŁA 80-4-57-54-00 	RĘKOJEŚĆ 80-28-5-1 	RĘKOJEŚĆ 80-28-5-1 
T25 12,14,16...90 90,95,100...115 L(mm) REF	80-X-01-04-L	TRZPIEŃ 	T25 14,16,18...60 65,70,75...95 L(mm) REF	80-X-01-17-L	TRZPIEŃ 	 S3.5 80-4-43-45-11	 S3.5 80-4-43-45-11
		 T25 80-43-268-T25.115			 T25 80-43-268-T25.115		