


**PŁYTKA -TPLO-KOMPRESYJNO-BLOKOWANA
 TYTAN**


L(mm)	Otwory	LEWA	PRAWA
25	6	80-1-25-241-25	80-1-25-254-25
30	6	80-1-25-241-30	80-1-25-254-30
35	6	80-1-25-241-35	80-1-25-254-35

Wkręty korowe samogwintujące Ø 1.5 Ø 1.1 T4	WIERTŁO Ø1.1 80-40-23-11.93.AO	
	PROWADNIK WIERTŁA 80-44-06-15	
	RĘKOJEŚĆ 80-101-10	
	TRZPIEŃ 80-43-270-T4.70	



Wkręty blokowane samogwintujące Ø 1.5 Ø 1.1 T4	WIERTŁO Ø1.1 80-40-23-11.93.AO	
	PROWADNIK WIERTŁA 80-44-06-15	
	RĘKOJEŚĆ 80-101-10	
	TRZPIEŃ 80-43-270-T4.70	

Wkręty blokowane zmiennokątowe Ø 1.5 Ø 1.1 T4	WIERTŁO Ø1.1 80-40-23-11.93.AO	
	PROWADNIK WIERTŁA 80-44-07-15	
	RĘKOJEŚĆ 80-101-10	
	TRZPIEŃ 80-43-270-T4.70	





Wkręty korowe samogwintujące Ø 2.0 Ø 1.5 T6	WIERTŁO Ø1.5 80-4-40-03-15	
	PROWADNIK WIERTŁA 80-4-57-191-00	
	RĘKOJEŚĆ 80-28-5-1	
	TRZPIEŃ 80-28-5-2 T6 80-43-268-T6.115	


Wkręty blokowane samogwintujące Ø 2.0 Ø 1.5 T6	WIERTŁO Ø1.5 80-4-40-03-15	
	PROWADNIK WIERTŁA 80-57-122-00	
	RĘKOJEŚĆ 80-28-5-1	
	TRZPIEŃ S1.5 80-28-5-2 T6 80-43-268-T6.115	

Wkręty blokowane zmiennokątowe Ø 2.0 Ø 1.5 T6	WIERTŁO Ø1.5 80-4-40-03-15	
	PROWADNIK WIERTŁA 80-4-57-212-00	
	RĘKOJEŚĆ 80-28-5-1	
	TRZPIEŃ 80-43-268-T6.115	

PŁYTKA -TPLO-KOMPRESYJNO-BLOKOWANA
 TYTAN


L(mm)	Otwory	LEWA	PRAWA
31,5	5	80-1-25-32-05	80-1-25-33-05
38	6	80-1-25-32-06/1	80-1-25-33-06/1
46	6	80-1-25-32-06	80-1-25-33-06

Wkręty korowe Ø 2.4 samogwintujące 		WIERTŁO Ø1.8 		Wkręty korowe Ø 2.7 samogwintujące 		WIERTŁO Ø2.0 	
T8	L(mm) REF 8,10,12...40 80-X-01-240-L	PROWADNIK WIERTŁA	80-4-57-192-00	T8	L(mm) REF 8,10,12...60 80-X-01-216-L	PROWADNIK WIERTŁA	80-4-57-192-00
	L(mm) REF 8,10,12...60 80-X-01-239-L	RĘKOJEŚĆ	80-28-5-1		L(mm) REF 8,10,12...60 80-X-01-48-L	RĘKOJEŚĆ	80-28-5-1
		TRZPIEŃ				TRZPIEŃ	
		○ S2.5	80-28-5-3			○ S2.5	80-28-5-3
		⊗ T8	80-43-268-T8.115			⊗ T8	80-43-268-T8.115
Wkręty blokowane Ø 2.4 samogwintujące 		WIERTŁO Ø1.8 		Wkręty blokowane Ø 2.7 samogwintujące 		WIERTŁO Ø2.0 	
T8	L(mm) REF 8,10,12...40 80-X-01-208-L	PROWADNIK WIERTŁA	80-4-57-49-00	T8	L(mm) REF 10,12,14...60 80-X-01-209-L	PROWADNIK WIERTŁA	80-4-57-74-00
	L(mm) REF 8,10,12...40 80-X-01-200-L	RĘKOJEŚĆ	80-28-5-1		L(mm) REF 8,10,12...60 80-X-01-201-L	RĘKOJEŚĆ	80-28-5-1
		TRZPIEŃ				TRZPIEŃ	
		○ S2.5	80-28-5-3			○ S2.5	80-28-5-3
		⊗ T8	80-43-268-T8.115			⊗ T8	80-43-268-T8.115
Wkręty blokowane zmiennokątowe Ø 2.4 		WIERTŁO Ø1.8 		Wkręty blokowane zmiennokątowe Ø 2.7 		WIERTŁO Ø2.0 	
T8	L(mm) REF 8,10,12...40 80-X-01-298-L	PROWADNIK WIERTŁA	80-4-57-171-01	T8	L(mm) REF 10,16,18...40 80-X-01-325-L	PROWADNIK WIERTŁA	80-4-57-171-02
		RĘKOJEŚĆ	80-28-5-1			RĘKOJEŚĆ	80-28-5-1
		TRZPIEŃ	80-43-268-T8.115			TRZPIEŃ	80-43-268-T8.115

PŁYTKA -TPLO-KOMPRESYJNO-BLOKOWANA
 TYTAN


L(mm)	Otwory	LEWA	PRAWA
60	6	80-1-25-30-06/1	80-1-25-31-06/1
72	6	80-1-25-30-06	80-1-25-31-06
88	8	80-1-25-30-08	80-1-25-31-08

Wkręty korowe Ø 3.5 samogwintujące Ø 2.5	
L(mm)	REF
T15 8,10,12...70	80-X-01-221-L
75,80,85...90	
L(mm)	REF
8,10,12...70	80-X-01-51-L
75,80,85...90	

WIERTŁO Ø2.5	
PROWADNIK WIERTŁA	80-4-57-191-00
RĘKOJEŚĆ	80-28-5-1
TRZPIEŃ	
S2.5	80-28-5-3
T15	80-43-268-T15.115


Wkręty blokowane Ø 3.5 samogwintujące Ø 2.7	
L(mm)	REF
T15 10,12,14...60	80-X-01-223-L
65,70,75...95	
L(mm)	REF
10,12,14...60	80-X-01-11-L
65,70,75...95	

WIERTŁO Ø2.7	
PROWADNIK WIERTŁA	80-4-57-74-00
RĘKOJEŚĆ	80-28-5-1
TRZPIEŃ	
S2.5	80-28-5-3
T15	80-43-268-T15.115

Wkręty blokowane zmiennokątowe Ø 3.5 Ø 2.7	
L(mm)	REF
T15 10,12,14...60	80-X-01-322-L
65,70,75...95	

WIERTŁO Ø2.7	
PROWADNIK WIERTŁA	80-4-57-171-03
RĘKOJEŚĆ	80-28-5-1
TRZPIEŃ	80-43-268-T15.115

PŁYTKA -TPLO-KOMPRESYJNO-BLOKOWANA
TYTAN




L(mm)	Otwory	LEWA	PRAWA
90	4	80-1-25-244-90	80-1-25-250-90
98	5	80-1-25-244-98	80-1-25-250-98


Wkręty korowe Ø 4.5 samogwintujące Ø 3.2 L(mm) REF		WIERTŁO Ø3.2 80-4-40-08-32		Wkręty blokowane Ø 5.0 samogwintujące Ø 3.7 L(mm) REF		WIERTŁO Ø3.7 80-4-40-06-37	
T25 12,14,16...90 80-X-01-218-L 90,95,100...115 L(mm) REF		PROWADNIK WIERTŁA 80-4-57-57-00		T25 14,16,18...60 80-X-01-222-L 65,70,75...95 L(mm) REF		PROWADNIK WIERTŁA 80-4-57-54-00	
80-X-01-04-L 90,95,100...115		RĘKOJEŚĆ 80-28-5-1		80-X-01-17-L 14,16,18...60 65,70,75...95		RĘKOJEŚĆ 80-28-5-1	
		TRZPIEŃ S3.5 80-4-43-45-11 T25 80-43-268-T25.115				TRZPIEŃ S3.5 80-4-43-45-11 T25 80-43-268-T25.115	

PŁYTKA-TPLO-PROSTA KOMPRESYJNA

TYTAN | STAL




L(mm)	t1	t2	Otworki	Nr Kat.
45	8	7	6	80-X-25-93-01
52	11	11	6	80-X-25-93-02



PŁYTKA-TPLO-PROSTA KOMPRESYJNA

TYTAN | STAL



L(mm)	t1	t2	Otworki	Nr Kat.
55	10	8	6	80-X-25-94-01
63	14	12	6	80-X-25-94-02
71	12	11	7	80-X-25-94-03
79	14	14	7	80-X-25-94-04

 Wkręty korowe Ø 2.4
 samogwintujące
 Ø 1.8

 L(mm) REF
T8 8,10,12...40 **80-X-01-240-L**

 L(mm) REF
 8,10,12...60 **80-X-01-239-L**

 WIERTŁO
 Ø1.8

 PROWADNIK
 WIERTŁA 80-4-57-192-00

RĘKOJEŚĆ 80-28-5-1

TRZPIEŃ

S2.5 80-28-5-3

T8 80-43-268-T8.115

 Wkręty korowe Ø 2.7
 samogwintujące
 Ø 2.0

 L(mm) REF
T8 8,10,12...60 **80-X-01-216-L**

 L(mm) REF
 8,10,12...60 **80-X-01-48-L**

 WIERTŁO
 Ø2.0

 PROWADNIK
 WIERTŁA 80-4-57-192-00

RĘKOJEŚĆ 80-28-5-1

TRZPIEŃ

S2.5 80-28-5-3

T8 80-43-268-T8.115

 Wkręty korowe Ø 3.5
 samogwintujące
 Ø 2.5

 L(mm) REF
T15 8,10,12...70 **80-X-01-221-L**

 L(mm) REF
 75,80,85...90

 L(mm) REF
 8,10,12...70 **80-X-01-51-L**

 L(mm) REF
 75,80,85...90

 WIERTŁO
 Ø2.5

 PROWADNIK
 WIERTŁA 80-4-57-191-00

RĘKOJEŚĆ 80-28-5-1

TRZPIEŃ

S2.5 80-28-5-3

T15 80-43-268-T15.115